



ULTRASONIC SCALER

Buyer's guide

We've collected all the ultrasonic scalers available in the UK so you can choose the best one for your practice.



Ultrasonic scalers use high frequency oscillations to clean teeth, gums, and implants. A unit costs on average £1,700, making it worth choosing one carefully. Key features to consider:

OSCILLATION TYPE

The method used to generate the oscillations, either piezoelectric or magnetostrictive. While little differs for the user, it is important for the buyer to note which type, as handpieces are not cross-compatible between the two.

PROCEDURES

Whether the manufacturer has designed the unit with specific procedures in mind, like endodontics.

CONTROLS

The type of controls, such as touchscreens, physical buttons, or rotating dials. Touchscreens are easier to clean and can be used while wearing surgical gloves. Physical buttons and dials are more tactile and generally cheaper.

FREQUENCY

The oscillation frequency of the unit. This is important for handpiece compatibility, and higher frequencies are quieter for the patient.

LED

Whether the handpiece has integrated LED lights to better illuminate the working area. Many units come in LED and non-LED versions.

AIR POLISHER

Some units have an integrated air polisher.

INTEGRATED WATER SUPPLY

Units with an integrated water supply are effectively standalone. Units without integrated water must be plumbed in.

NB 20240808-2

	TYPE	PROCEDURES	CONTROLS	FREQUENCY (KHZ)	LED	AIR POLISH-ER	INTEGRATED WATER SUPPLY	PRICE INC VAT
--	------	------------	----------	-----------------	-----	---------------	-------------------------	---------------

DTE - China

UK distributors: globaldental.co.uk, nextdental.com, dentalwarehouse.co.uk, ukdentalsupplies.com, qudent.co.uk, mcdental.co.uk, dentaledgeuk.co.uk, presidential.co.uk, toothsaver.co.uk, vsdent.com, medident-uk.com, f2medicalsupplies.com, dentaldeal.co.uk

D5	Piezo	Periodontics, endodontics	Buttons	28±3	Yes (D5 LED variant)	No	No	£316
D600	Piezo	Periodontics, endodontics	Touch panel	28±3	Yes (D600 LED variant)	No	700ml	£654
D7	Piezo	Periodontics, endodontics	Buttons	28±3	Yes (D7 LED variant)	No	500ml	£349
S6	Piezo	Periodontics, endodontics	Touch panel	28±3	Yes (S7 LED variant)	No	No	£352

Woodpecker - China

UK distributors: qudent.co.uk, medimatch.co.uk, f2medicalsupplies.com, nextdental.com, mcdental.co.uk, dentaledgeuk.co.uk, dentalwarehouse.co.uk, ukdentalsupplies.com, trycare.co.uk, dentaldeal.co.uk, dentalsky.com, toothsaver.co.uk, southwestmedical.co.uk, presidential.co.uk

AP-B	Piezo	Periodontics, polishing	Buttons	30±5	Yes	Yes	700ml & 1400ml	£2,280
PT-A	Piezo	Periodontics, endodontics, implant maintenance, polishing	Touchscreen	30±5	No	Yes	600ml & 1400ml	£3,200
PT-B	Piezo	Periodontics, endodontics, implant maintenance, polishing	Touchscreen	30±5	No	Yes	600ml & 1400ml	£5,160
U6	Piezo	Periodontics, endodontics	Touch panel	28±3	Yes (U6 LED variant)	No	No	£265
U600	Piezo	Periodontics, endodontics	Touch panel	28±3	Yes (U600 LED variant)	No	700ml	£475
UDS-A	Piezo	Periodontics, endodontics	Buttons	28±3	Yes (UDS-A LED variant)	No	No	£221
UDS-E	Piezo	Periodontics, endodontics	Buttons	28±3	Yes (UDS-E LED variant)	No	350ml	£347
UDS-J	Piezo		Dial	28±3	No	No	No	£270
UDS-K	Piezo		Dial	28±3	Yes (UDS-K LED variant)	No	No	£204
UDS-L	Piezo	Periodontics, endodontics	Buttons	28±3	No	No	900ml	£286
UDS-P LED	Piezo	Periodontics, endodontics	Buttons	28±3	Yes	No	No	£270

	TYPE	PROCEDURES	CONTROLS	FREQUENCY (KHZ)	LED	AIR POLISH-ER	INTEGRATED WATER SUPPLY	PRICE INC VAT
--	------	------------	----------	-----------------	-----	---------------	-------------------------	---------------

Bonart - Taiwan

UK distributors: dentaledgeuk.co.uk, mcdental.co.uk, nextdental.com, f2medicalsupplies.com, ukdentalsupplies.com

Marquee ART-M3II	Magneto	Periodontics, endodontics	Dial	25 or 30	No	No	No	£510
Piper ART-P6	Piezo		Dial	26-32	No	No	No	£510

EMS - Switzerland

UK distributors: mcdental.co.uk, trycare.co.uk, haguedental.com

Piezon 250	Piezo	Periodontics	Dial	24-32	Yes	No	350ml or 500ml	£2,038
Piezon Master 700	Piezo	Endodontics	Touch panel	24-32	Yes	No	350ml & 500ml	£4,314
Airflow prophylaxis master	Piezo	Periodontics, endodontics, polishing	Touch panel	24-32	No	Yes	800ml	£9,733

W&H - Austria

UK distributors: mcdental.co.uk, swallowdental.co.uk

Tigon+	Piezo	Periodontics, endodontics	Buttons	27-32	Yes	No	400ml	£1,426
Pyon 2	Piezo	Periodontics, endodontics	Dial	27-32	Optional	No	No	£764
Proxeo Ultra	Piezo	Periodontics, endodontics	Dial	22-35	Yes	No	750ml (PB-520 & PB-530 only)	£1,084

NSK-Nakanishi - Japan

UK distributors: ukdentistry.co.uk, dentaledgeuk.co.uk, swallowdental.co.uk, ukdentalsupplies.com

Varios 970	Piezo	Periodontics, endodontics	Touch panel	28-32	Yes (the Optic variant)	No	2 x 400ml	£1,495
Varios Combi Pro	Piezo	Polishing	Touch panel	28-32	Yes (the Optic variant)	Yes	2 x 400ml	£4,756
Varios 570	Piezo		Dial	28-32	Yes	No	No	£1,195
Varios 370	Piezo		Dial	28-32	Yes (the Optic variant)	No	No	£725

	TYPE	PROCEDURES	CONTROLS	FREQUENCY (KHZ)	LED	AIR POLISH-ER	INTEGRATED WATER SUPPLY	PRICE INC VAT
--	------	------------	----------	-----------------	-----	---------------	-------------------------	---------------

Mectron - Italy

UK distributors: mcdental.co.uk

Micropiezo S	Piezo	Periodontics	Buttons	24-36	Yes	No	No	£625
--------------	-------	--------------	---------	-------	-----	----	----	------

Acteon/Satelec - France

UK distributors: qudent.co.uk, henryschein.co.uk, cadentltd.co.uk, presidential.co.uk, globaldental.co.uk, dental-equipment.co.uk, ssdgroup.co.uk, ukdentalsupplies.com

Newtron booster	Piezo		Dial	>28	No	No	No	£1,333
Newtron P5 XS BLED	Piezo		Dial	>28	Yes	No	300ml	£2,048
Newtron P5 BLED	Piezo		Dial	>28	Yes	No	No	£1,323

VRN - China

UK distributors: dentaledgeuk.co.uk

A8	Piezo	Periodontics, endodontics	Buttons	25-31	Yes	No	500ml	£432
DA-20	Piezo	Periodontics, endodontics	Touch panel	18-45	Yes	No	800ml	£660
Q6	Piezo	Periodontics, endodontics, polishing	Touch panel	28±3	No	Yes	2 x 600ml	£2,754

	TYPE	PROCEDURES	CONTROLS	FREQUENCY (KHZ)	LED	AIR POLISH-ER	INTEGRATED WATER SUPPLY	PRICE INC VAT
--	------	------------	----------	-----------------	-----	---------------	-------------------------	---------------

Dentsply Sirona - USA

UK distributors: dentalsky.com, dkggroup.com, kentexpress.co.uk, bfmulholland.com, trycare.co.uk, nextdental.com

Cavitron 300	Magneto		Touchscreen	30	No	No	No	£2,994
Cavitron Plus	Magneto		Buttons	30	No	No	No	£3,054
Cavitron Select SPS	Magneto		Dial	30	No	No	Optional	£3,054
Cavitron Touch	Magneto		Touchscreen	30	No	No	No	£3,766

Rolence - Taiwan

UK distributors: ukdentalsupplies.com, dentalwarehouse.co.uk

Elitedent MS-1 Dual	Magneto	Periodontics, endodontics	Dial	25-30	No	No	2 x 500ml	£719
---------------------	---------	---------------------------	------	-------	----	----	-----------	------

R&S - France

UK distributors: dentalsky.com, clinipline.co.uk

Tri Scaler Aqua	Piezo	Periodontics, endodontics	Buttons		No	No	300ml or 500ml	£545
Tri Scaler Compact	Piezo	Periodontics, endodontics	Buttons		No	No	No	£330

3TECH - France

UK distributors: dentalsky.com, clinipline.co.uk

Easy scaler L-28	Piezo	Periodontics, endodontics	Touch panel	28±3	Yes	No	600ml	£432
------------------	-------	---------------------------	-------------	------	-----	----	-------	------

B.A. International - UK

UK distributors: dentalsky.com, clinipline.co.uk

Ulticlean UC500L	Piezo	Periodontics, endodontics, implant maintenance, polishing	Touchscreen	30±5	Yes	Yes	1500ml & 600ml	£5,519
------------------	-------	---	-------------	------	-----	-----	----------------	--------

Questions? Comments?

We'd love to hear from you! Contact us at find@eltoco.com

Stop searching

START FINDING

