



X-RAY PLATE SCANNER

Buyer's guide

Choose the best x-plate plate scanner for your practice with our handy guide to all the models available in the UK today.



Plate scanners digitise oral x-rays for later storage and analysis. They range widely in price, depending on brand and functionality. Key features to consider:

VIEWSCREEN

An integrated viewscreen allows x-rays to be viewed without opening them on a computer.

STANDALONE MODE

Some scanners also store images on the device, allowing the practitioner to work entirely independently of a computer.

MULTIPLE PLATES

Scanners which can read multiple plates simultaneously are larger but more efficient.

CONTROL TYPE

The type of controls on the scanner unit.

PLATE SIZES

The range of x-ray plate sizes which the scanner can handle. Many scanners cannot read size 4 plates.

CONNECTIVITY

How the scanner can connect to a computer to upload and store the digitised images. Ethernet connections may require running a cable to your router so if this a problem WIFI or USB is preferable.

RESOLUTION

The resolution of the scanner, measured in line pairs per millimetre. Higher is better. Actual resolutions are lower than this but not comparable across manufacturers due to differences in how this is calculated.

NB 20241101-1

	VIEWSCREEN	STANDALONE MODE	MULTIPLE PLATES	CONTROLS	PLATE SIZES	CONNECTIVITY	RESOLUTION (LP/MM)	SIZE W X H X D (CM)	PRICE INC VAT
--	------------	-----------------	-----------------	----------	-------------	--------------	--------------------	---------------------	---------------

Acteon/Satelec - France

UK distributors: ukdentalsupplies.com, 360v.co.uk, ssdgroup.co.uk, dentalsky.com, wysdom.co.uk

PSPIX2	Yes	No	No	Touchscreen	0-4	Ethernet	20	20x15x19	£4,794
--------	-----	----	----	-------------	-----	----------	----	----------	--------

Carestream - USA

UK distributors: 360v.co.uk, hulbertimaging.co.uk, thedentalimagingcompany.co.uk, dmi.co.uk

CS 7200	No	No	No	Buttons	0-2	Ethernet & USB	19	13x27x30	£5,039
CS 7600	No	No	No	Touchscreen	0-4	Ethernet	18	23.7x26.7x26	£5,999

Dexis - USA

UK distributors: dentalimaging.co.uk, angliandental.co.uk, dental dynamix.co.uk, midental.co.uk

Scan eXam	No	No	No	Buttons	0-3	Ethernet	17	21.5x19.6x38.2	
Scan eXam One	No	No	No	Buttons	0-3	Ethernet	17	23.7x26.7x26	£4,620

Dürr Dental - Germany

UK distributors: dentalimaging.co.uk, 360v.co.uk, probomedical.co.uk, dental edgeuk.co.uk, agdental.co.uk, haguedental.com, bfmulholland.com, ssdgroup.co.uk, dentalsky.com, dental-equipment.co.uk, dbdental.co.uk, woodlanedental.co.uk, sgdentalsupplies.co.uk, mckillopdental.co.uk, dentalengineers.com, avensysmedical.co.uk, bowendentaluk.com, wysdom.co.uk

VistaScan Ultra View	Yes	Yes	Yes	Touchscreen	0-4	WIFI & Ethernet	40	25.1x39.1x29.2	£10,930
VistaScan Mini View	Yes	Yes	No	Touchscreen	0-4	WIFI & Ethernet	40	22.6x27.5x24.3	£8,159
VistaScan Mini View 2.0	Yes	Yes	No	Touchscreen	0-4	WIFI & Ethernet	40	21.1x27.3x25.8	£10,199
VistaScan Mini Easy	No	No	No	Buttons	0, 2	Ethernet & USB	40	23.4x22.6x24.3	£8,184
VistaScan Mini Easy 2.0	No	No	No	Touchscreen	0, 2 (1, 3, 4 optional)	Ethernet	40	21.1x24.9x25.8	£6,044
VistaScan Combi View	Yes	Yes	Yes	Touchscreen	0-4	WIFI & Ethernet	40	38x72x38	£10,192
VistaScan Nano Easy	No	No	No	Buttons	2 (0, 1 optional)	Ethernet	17	16.7x23.1x21.6	

	VIEWSCREEN	STANDALONE MODE	MULTIPLE PLATES	CONTROLS	PLATE SIZES	CONNECTIVITY	RESOLUTION (LP/MM)	SIZE W X H X D (CM)	PRICE INC VAT
--	------------	-----------------	-----------------	----------	-------------	--------------	--------------------	---------------------	---------------

NewTom - Italy

UK distributors: hsdequipment.co.uk

X-PSP	No	No	No	None - fully automated	0-3	Ethernet	17	16.3x22.4x29	
-------	----	----	----	------------------------	-----	----------	----	--------------	--

Planmeca - Finland

UK distributors: dental dynamix.co.uk, dental edgeuk.co.uk, aihealthcare.uk, currandental.com, wysdom.co.uk, haguedental.com, protechdental.co.uk, deca-dent.co.uk

ProScanner 2.0	No	No	No	Buttons	0-3	Ethernet	17	16.7x23.1x21.6	
----------------	----	----	----	---------	-----	----------	----	----------------	--

Woodpecker - China

UK distributors: f2medicalsupplies.com, medimatch.co.uk, vsdent.com, toothsaver.co.uk

I Scan	Yes	Yes	No	Touchscreen	0-3	WIFI & USB	17	13.8x24.6x29	£2,999
Free Scan	No	No	No	Touch Button	0-3	USB	17	16.5x21x28	£2,878

Questions? Comments?

We'd love to hear from you! Contact us at find@eltoco.com

Stop searching
START FINDING

